



104 Shorting Road
Toronto, Ontario
M1S 3S4 Canada
Telephone 416.299.4003
Toll Free 800.463.9613
Facsimile 416.299.4912

October 2, 2012

To whom it may concern

Statement regarding Hyperbaric Oxygen Therapy

The following Derma Products sold in the USA, may be used in conjunction with Hyperbaric Oxygen Therapy.

Medihoney® Calcium Alginate Dressing

Medihoney® Honeycolloid Dressing

Medihoney® Paste

Medihoney® Gel

Medihoney® HCS (Hydrogel Colloidal Sheet) Dressing

Xtrasorb™ Super Absorbent Dressing

Xtrasorb™ Foam Super Absorbent Dressing

Xtrasorb™ HCS (Hydrogel Colloidal Sheet) Dressing

TCC-EZ® (must use in HBO treatment, only after 24 hrs of curing of the cast)

MEDE-KAST™ (must use in HBO treatment, only after 24 hrs of curing of the cast)

MEDE-KAST™ ULTRA (must use in HBO treatment, only after 24 hrs of curing of the cast)

Based on the flash point studies and on a risk analysis performed of each product and on feedback that we currently receive from our customer complaint system, we conclude that there are no additional risks when using the above products in Hyperbaric chambers.

This declaration is made to the best of our knowledge and may be reconsidered if new facts turn up.

Best Regards,

Sharmini Atheray

VP, Quality & Regulatory Affairs